Description

Tetric® N-Ceram Bulk Fill is a state-of-the-art, light-curing, radiopaque hybrid composite for direct restorations in posterior teeth. Tetric N-Ceram Bulk Fill cures with light in the wavelength range of 400–500 nm (blue light) and can be applied in layers of up to 4 mm.

Shades

Tetric N-Ceram Bulk Fill is available in the following 3 shades:

- Universal shades: IVA, IVB, IVW

Composition

The monomer matrix is composed of dimethacrylates (19–21% weight). The total content of inorganic fillers is 75–77% weight or 53–55% volume. The fillers consist of barium glass, prepolymer, ytterbium trifluoride and mixed oxide. Additives, catalysts, stabilizers and pigments are additional contents (< 1.0% weight). The particle size of the inorganic fillers is between 0.04 and 3 µm. The mean particle size is 0.6 µm.

Indication

- Restoration of deciduous teeth
- Restorations in the posterior region (Classes I and II, including the replacement of individual cusps)
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Reconstructive build-up
- Extended fissure sealing in molars and premolars

Contraindication

Placement of Tetric N-Ceram Bulk Fill restorations is contraindicated

- if a dry working field cannot be established or the prescribed application technique cannot be applied;
- if the patient is known to be allergic to any of the ingredients in Tetric N-Ceram Bulk Fill.

Side effects

In individual cases, components of Tetric N-Ceram Bulk Fill may lead to sensitization. The product should not be used in such cases. To avoid possible irritation of the pulp, areas close to the pulp should be protected with a suitable pulp/dentin protector (selectively apply a calcium hydroxide-based preparation in areas close to the pulp and cover with suitable cavity liner).

Interactions

Phenolic substances such as eugenol/clove oil inhibit the polymerization of materials based on methacrylate. Consequently, the application of such materials together with Tetric N-Ceram Bulk Fill must be avoided. Discolouration may occur in combination with cationic mouthwashes, plaque disclosing agents and chlorhexidine.

Usage**1. Shade determination**

Clean the teeth prior to shade determination. The shade is selected with the tooth still moist. To facilitate the procedure, the Tetric N-Ceram Bulk Fill shade guide can be used.

2. Isolation

Appropriate isolation is required, preferably with a rubber dam (e.g. OptraDam® Plus).

3. Cavity preparation

The cavity is prepared according to the principles of the adhesive technique, i.e., by preserving as much of the tooth structure as possible. Do not prepare sharp internal edges or angles or additional undercuts in caries-free areas. The dimensions of the cavity are generally dictated by the extent of the caries or the size of the previous restorations. Bevel enamel edges of anterior teeth. In the posterior region, only the sharp enamel edges should be slightly broken or rounded (finishing diamonds, 25–40 µm). Caries-free cervical defects are not prepared, only cleaned with pumice or other suitable cleaning pastes with the help of rubber polishers or rotary brushes. Subsequently, rinse the cavity with water to remove all the residue and dry it with water- and oil-free air.

4. Pulp protection / Base

Do not apply a base if an enamel-dentin bonding agent is used. In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide-based preparation (e.g. ApexCal®) and subsequently covered with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner). Do not cover the remaining cavity walls so that they can be used to generate a bond with an enamel-dentin adhesive.

5. Apply matrix / interdental wedge

Use a wrap around matrix for cavities affecting the proximal area or a sectional matrix and wedge it.

6. Conditioning / Application of the bonding agent

Condition and apply the bonding agent according to the instructions for use of the product in use. Ivoclar Vivadent recommends using Tetric® N-Bond (with phosphoric acid etching, e.g. N-Etch) or the self-etching adhesive Tetric® N-Bond Self-Etch.

7. Application of Tetric N-Ceram Bulk Fill

- In order to achieve optimum results, Tetric N-Ceram Bulk Fill should be applied in increments of max. 4 mm and adapted to the cavity walls with a suitable instrument (e.g. OptraSculpt®). Use a special contact point instrument (e.g. the bifurcated OptraContact®) for large cavities where contacts are difficult to create.
- Sufficient exposure to the curing light prevents incomplete polymerization. For the recommendations regarding exposure time and light intensity see Table 1.
- When using a metal matrix, additionally polymerize the composite material from the buccal or the lingual/palatal aspect after removing the matrix, if no Bluephase® polymerization light is used or the light probe cannot be ideally positioned, e.g. distant to composite or diverging scattering angle.

Tetric® N-Ceram Bulk Fill

EN Instructions for Use

Light-curing resin-based dental restorative material

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- In many cases, a flowable composite is used today as an initial layer to create an even cavity floor and to facilitate the adaptation of the subsequently used restorative material. A flowable composite (e.g. Tetric® N-Flow) can be applied in a thin initial layer. This optional layer must be cured separately (please refer to the respective Instructions for Use).

8. Finishing / Checking the occlusion / Polishing

After polymerization, remove excess material with suitable finishers or fine diamonds. Remove proximal excess with diamond or tungsten carbide finishers or finishing strips or flexible finishing discs. Check the occlusion and articulation and grind in the necessary adjustments to prevent premature contacts or undesired occlusal paths on the surface of the restoration. Use silicone polishers (e.g. OptraPol® Next Generation) as well as polishing discs and polishing strips to polish the restoration to a high gloss.

Additional notes

1. Tetric N-Ceram Bulk Fill can be used in combination with Tetric N-Ceram and Tetric N-Flow.
2. In the case of repairs, additional Tetric N-Ceram Bulk Fill can be directly applied to polymerized material. If the Tetric N-Ceram Bulk Fill restoration has already been polished, it has to be roughened and wetted with Heliobond before a new layer of Tetric N-Ceram Bulk Fill can be applied.
3. Tetric N-Ceram Bulk Fill should have ambient temperature when applied. Cool temperatures render the material difficult to extrude.
4.  For single use only. If Tetric N-Ceram Bulk Fill is applied from the Cavifil directly in the mouth of the patient, the Cavifil must not be used for more than one patient due to hygienic reasons (prevention of cross-contamination between patients).
5. Do not use oxidizing disinfectants to disinfect syringes and Cavifils.
6. The recommended increment thickness is based on hardness profile measurement.

Warning

Avoid contact of unpolymerized Tetric N-Ceram Bulk Fill with the skin/mucous membrane and eyes. Unpolymerized Tetric N-Ceram Bulk Fill may have a slight irritating effect and may lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Shelf life and storage

- Do not use Tetric N-Ceram Bulk Fill after the expiry date.
- Storage temperature 2 – 28 °C (36 – 82 °F).
- Close syringes/Cavifils immediately after usage. Exposure to light causes premature polymerization.
- Expiry date: see note on Cavifils, syringes and packages.

Table 1

UNIT	Bluephase C8	Bluephase	Bluephase 20i	Bluephase Style
PROGRAM				
Turbo	–	–	5 s	–
High Power	15 s	10 s	10 s	10 s
Soft Start	20 s	15 s	15 s	–

Light Intensity	Exposure Time
≥ 500 mW/cm ²	20 s
≥ 1'000 mW/cm ²	10 s

Keep material out of children's reach.

For use in dentistry only.

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Pte. Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.